

M1030
-06-36


Novel Biphasic Lipolysis Testing of Nilotinib Lipid Based Suspensions to Predict *In Vivo* Performance.

Patrick J. O'Dwyer^{1,2}; Niklas J. Koehl³, Karl J. Box¹; Brendan T. Griffin³, René Holm⁴, Martin Kuentz⁵, Christos Reppas²

¹ Research and Development, Pion Inc. (UK) Ltd, Forest Row RH18 5DW, UK; ² Department of Pharmacy, School of Health Sciences, National and Kapodistrian University of Athens, Zografou, Greece; ³ School of Pharmacy, University College Cork, College Road, Cork, Ireland; ⁴ Drug Product Development, Janssen Research and Development Johnson & Johnson Turnhoutseweg 30, 2340 Beerse, Belgium; ⁵ Institute of Pharma Technology University of Applied Sciences and Arts Northwestern Switzerland 4132 Muttenz, Switzerland.



Advancing Pharmaceutical Sciences,
Careers, and Community



CONTACT INFORMATION: podwyer@pion-inc.com

PURPOSE

- Due to the low aqueous solubility of many new chemical entities, lipid based formulations (LBF) have generated significant interest as a way to increase oral bioavailability.
- However, suitable *in vitro* lipolysis testing remains challenging, with many methods producing results which are quite distinct from *in vivo* performance [1].

OBJECTIVE(S)

- The objective of this study was to develop a new and quick two stage *in vitro* lipolysis protocol to incorporate the absorption of drug in the small intestine, using a decanol layer.
- Nilotinib, a weakly basic drug, was selected as an appropriate model 'brick-dust' drug. Results were compared to a standard *in vitro* lipolysis experiment (pH-stat method) and the results from an *in vivo* rat study [2].

METHOD(S)

- Four lipid suspensions were tested containing crystalline nilotinib. The biphasic lipolysis was carried out using the inForm (Pion Inc.) instrument with the experimental setup shown in Figure 1.
- At the end of the experiment in the intestinal sector, the pH was back to pH 9 to determine the release of non-ionized free fatty acids.
- Drug was quantified *in situ* in real time by fibre optic UV dip probes.
- A higher stirring rate (300rpm) was applied during the gastric sector and during the dispersion period with the pancreatic enzymes (compared to a stirring rate of 100rpm during the intestinal sector). This was done to ensure the formulations were well dispersed prior to the addition of the decanol layer and prevent 'non-physiological' mixing of the LBFs with the decanol layer.

RESULT(S)

- Data from the standard *in vitro* method (pH-stat method) indicated a greater release from the Capmul MCM formulation than the Peceol formulation. This was in contrast to the results from the clinical study and biphasic lipolysis experiment (Table 1).
- The standard *in vitro* method showed a significant difference in the release of the Captex and Olive Oil formulations, unlike the *in vivo* and biphasic lipolysis results (Table 1).
- In addition, the low concentration of drug in solution after 60 minutes (<1% of the dose for all formulations) raises concerns about the pH-stat method.

Table 1. Overview of the *in vivo* and *in vitro* results for the Nilotinib formulation tested. Each data point represents the mean \pm SD.

Formulation	Sprague Dawley rat AUC _{0-inf,h} [ng*h/mL] (n=5) [2]	Standard Lipolysis Method: Percentage Released in Aqueous Phase after 60 minutes Digestion (n=3) [2]	Biphasic Lipolysis Method: Percentage Released in Decanol Layer after 60 minutes Digestion (n=3)
Peceol	13,103 \pm 2557	0.51 \pm 0.06	36.21 \pm 1.49
Capmul MCM	11,210 \pm 5476	0.99 \pm 0.06	29.53 \pm 1.01
Captex 1000	5,168 \pm 2197	0.82 \pm 0.09	20.09 \pm 2.62
Olive Oil	3,548 \pm 2711	0.22 \pm 0.03	24.18 \pm 1.13

CONCLUSION(S)

Results from the pH-stat method did not provide any reasonable insight into the respective *in vivo* performance.

Compared to the pH-stat method, this novel biphasic lipolysis method offered:

- an improved prediction of *in vivo* performance of nilotinib LBFs
- avoided the requirement for centrifugation and HPLC analysis
- facilitated a more rapid *in vitro* screening process of novel LBFs.

Further assessment of this method will be carried out using different types of LBFs with clinical data from other species.

REFERENCE

1. Griffin BT, Kuentz M, Vertzoni M, Kostewicz ES, Fei Y, Faisal W, et al. Comparison of *in vitro* tests at various levels of complexity for the prediction of *in vivo* performance of lipid-based formulations: Case studies with fenofibrate. Eur J Pharm Biopharm [Internet]. 2014;86(3):427–37. Available from: <http://dx.doi.org/10.1016/j.ejpb.2013.10.016>
2. Koehl NJ, Holm R, Kuentz M, Griffin BT. New Insights into Using Lipid Based Suspensions for 'Brick Dust' Molecules : Case Study of Nilotinib. Pharm Res. 2019;36(56).

FUNDING

This project has received funding from the European Union's Horizon 2020 Research and Innovation Programme under grant agreement No 674909 (PEARRL)




Fig 1. Schematic of the Lipolysis Biphasic Dissolution setup.